



Optimam® LP/T



Optimam® MP/T



Optimam® HP/T



Optimam® XP/T

Optimam® *SublimeLine*®

POLYtxt®

Silicone-gel-filled mammary implants, with textured surface, oblong oval base, anatomical profile, maximum shape control

Implants made by POLYTECH – QUALITY made in Germany



POLYTECH
health & aesthetics 

Optimam® – POLYtxt®													
		LP Low Projection			MP Moderate Projection			HP High Projection			XP Extra-High Projection		
Base (mm)	Order n°	Projection (mm)	Vol. (ml)	Order n°	Projection (mm)	Vol. (ml)	Order n°	Projection (mm)	Vol. (ml)	Order n°	Projection (mm)	Vol. (ml)	
A B		C D			C D			C D			C D		
73 84	20774-060 [*]	23 35	60	20775-075 [*]	28 39	75	20776-095 [*]	34 44	95	20777-110 [*]	40 48	110	
78 90	20774-070 ¹	24 37	70	20775-090 ¹	30 41	90	20776-110 ¹	36 46	110	20777-130 ¹	43 51	130	
83 96	20774-085 [*]	24 39	85	20775-105 [*]	31 44	105	20776-130 [*]	37 49	130	20777-160 [*]	45 54	160	
88 102	20774-100 ¹	25 41	100	20775-125 ¹	32 46	125	20776-155 ¹	39 52	155	20777-185 ¹	47 57	185	
93 107	20774-110 [*]	26 43	110	20775-150 [*]	33 48	150	20776-180 [*]	40 54	180	20777-220 [*]	49 60	220	
98 113	20774-130 ¹	27 45	130	20775-170 ¹	35 51	170	20776-215 ¹	42 57	215	20777-255 ¹	51 63	255	
103 119	20774-150 [*]	28 47	150	20775-200 [*]	36 53	200	20776-250 [*]	44 60	250	20777-290 [*]	53 66	290	
108 125	20774-170 ¹	29 49	170	20775-225	37 56	225	20776-280	45 62	280	20777-335	55 69	335	
113 131	20774-195 [*]	29 51	195	20775-260 [*]	38 58	260	20776-320 [*]	47 65	320	20777-385 [*]	57 72	385	
118 136	20774-220 ¹	30 54	220	20775-290	40 60	290	20776-365	49 68	365	20777-440	59 75	440	
123 142	20774-240 [*]	31 56	240	20775-330 [*]	41 63	330	20776-415 [*]	51 70	415	20777-495 [*]	61 78	495	
128 148	20774-275 ¹	32 58	275	20775-370 ¹	43 65	370	20776-460 ¹	53 73	460	20777-555 ¹	63 81	555	
133 154	20774-305 [*]	33 60	305	20775-405 [*]	44 67	405	20776-515 [*]	55 76	515	20777-620 [*]	65 84	620	
138 160	20774-335 ¹	34 62	335	20775-455 ¹	45 70	455	20776-570 ¹	57 78	570	20777-695 ¹	67 87	695	
143 165	20774-375 [*]	35 64	375	20775-505 [*]	47 72	505	20776-640 [*]	59 81	640	20777-770 [*]	69 90	770	
148 171	20774-415 ¹	35 66	415	20775-555 ¹	48 74	555	20776-705 ¹	61 84	705	20777-855 ¹	71 93	855	
153 177	20774-450 [*]	36 68	450	20775-610 [*]	49 77	610	20776-775 [*]	62 86	775	20777-940 [*]	73 96	940	
158 183	20774-495 ¹	37 70	495	20775-670 ¹	51 79	670	20776-855 ¹	64 89	855	20777-1035 ¹	75 99	1035	

Approximate measurements. Sizes highlighted in color and marked (1) will be produced on demand; supply ≈ 6 weeks after order. Silicone gel-filled, smooth sizers (sterile, re-sterilisable) are available for products marked with an asterisk (*). Further sizes will be produced according to our terms for custom-made products.

Product Description

The silicone-gel-filled Optimam® *SublimeLine*® implants with POLYtxt® surface have an oblong oval base and an anatomical projection with the highest point in the lower third. The filling level of the *SublimeLine*® implants is maximized to provide more volume.

The consistency of these high-quality mammary implants is very close to the haptics and motional properties of the mammae. The implants are available in four projections: Low Profile (LP), Moderate Profile (MP), High Profile (HP), and Extra High Profile (XP).

The special curve measurement D, which describes the distance between the highest point of projection and the caudal implant edge, establishes a reference parameter for the correct placement of the areola in relation to the inframammary fold.

Information on Material

POLYTECH Health & Aesthetics implants are filled with the highly cross-linked EasyFit™ gel, which stands out due to its natural consistency. Because of the gel's "memory effect", it always returns to its original shape after moderate mechanical deformation.

The chemically and mechanically resistant, elastic shell of the implants consists of several layers of a silicone elastomer. In order to maintain a high level of safety, the

shell of the silicone gel-filled implants is equipped with an additional diffusion barrier layer.

The surface of the implant shell is POLYtxt® textured. This surface is known for its reliability and has been proven to contribute to a reduction of capsular contracture rates¹.

Production of the implants is based on a quality management according to ISO 9001 and ISO 13485. All materials used are of implantable medical quality.

¹: textured implants contracture rates: 8–15% (Iwuagwu + Frame, 1997), which is 15% less than with smooth implants (Handel, 2006)

Indications

- Breast reconstruction after mastectomy or ablatio mammae
- Contour augmentation and correction
- Adjustment of asymmetries, aplasia or hypoplasia
- Combined breast and thoracic deformities
- Ptosis or atrophy
- Implant exchange

Important Information

The implants are delivered individually packed and sterile.

Attention: Please note the detailed instructions for use included with each product.



POLYTECH Health & Aesthetics GmbH

Altheimer Str. 32 • 64807 Dieburg • Germany

+49 (0) 60 71 98 63 0 • +49 (0) 60 71 98 63 30

info@polytechhealth.com • www.polytech-health-aesthetics.com