



Même® LP/M



Même® MP/M



Même® HP/M



Même® XP/M

Même® *SublimeLine*®

Microthane®

Silicone-gel-filled mammary implants, with Microthane® surface, round base, central profile, maximum shape control

Implants made by POLYTECH – QUALITY made in Germany



Même® Microthane®													
		LP Low Projection			MP Moderate Projection			HP High Projection			XP Extra-High Projection		
Base Width (mm)	Order n°	Projection (mm)		Vol. (ml)	Order n°	Projection (mm)		Vol. (ml)	Order n°	Projection (mm)		Vol. (ml)	
A	C	B	D		B	D			B	D	B	D	
73	73	30724-055*	25 46	55	30725-070*	30 51	70	30726-090*	36 55	90	30727-105*	42 59	105
78	78	30724-065 ¹	26 49	65	30725-090 ¹	32 54	90	30726-110 ¹	38 58	110	30727-125	45 63	125
83	83	30724-085*	26 51	85	30725-105*	33 57	105	30726-125*	39 62	125	30727-155*	47 66	155
88	88	30724-095 ¹	27 54	95	30725-120	34 60	120	30726-155	41 65	155	30727-180	49 70	180
93	93	30724-110*	28 57	110	30725-145*	35 63	145	30726-180*	42 68	180	30727-215*	51 73	215
98	98	30724-125	29 60	125	30725-165	37 66	165	30726-210	44 71	210	30727-250	53 77	250
103	103	30724-150*	30 62	150	30725-195*	38 69	195	30726-235*	46 75	235	30727-285*	55 81	285
108	108	30724-165	31 65	165	30725-220	39 72	220	30726-275	47 78	275	30727-330	57 84	330
113	113	30724-185*	31 68	185	30725-255*	40 75	255	30726-315*	49 81	315	30727-380*	59 88	380
118	118	30724-215	32 71	215	30725-285	42 78	285	30726-360	51 84	360	30727-430	61 91	430
123	123	30724-235*	33 73	235	30725-320*	43 81	320	30726-400*	53 87	400	30727-485*	63 95	485
128	128	30724-270	34 76	270	30725-360	45 84	360	30726-450	55 91	450	30727-545	65 98	545
133	133	30724-295*	35 79	295	30725-400*	46 87	400	30726-505*	57 94	505	30727-610*	67 102	610
138	138	30724-330	36 82	330	30725-445	47 90	445	30726-560	59 97	560	30727-675	69 105	675
143	143	30724-365*	37 85	365	30725-495*	49 93	495	30726-620*	61 100	620	30727-755*	71 109	755
148	148	30724-400 ¹	37 87	400	30725-545 ¹	50 96	545	30726-690 ¹	63 104	690	30727-830 ¹	73 112	830
153	153	30724-440*	38 90	440	30725-600*	51 99	600	30726-760*	64 107	760	30727-920*	75 116	920
158	158	30724-485 ¹	39 93	485	30725-660 ¹	53 102	660	30726-835 ¹	66 110	835	30727-1010 ¹	77 119	1010

Approximate measurements. Sizes marked 1 and set in grey as well as further sizes will be produced on demand. Silicone gel-filled, smooth sizers (non sterile, re-sterilisable) are available for products marked with an asterisk *.

Product Description

The silicone-gel-filled Même® *SublimeLine*® implants with Microthane® surface have a round base and a central projection. The filling level of the *SublimeLine*® implants is maximized to provide more volume.

The consistency of these high-quality mammary implants is very close to the haptics and motional properties of the mammae. The implants are available in four projections: Low Profile (LP), Moderate Profile (MP), High Profile (HP), and Extra High Profile (XP).

The special curve measurement D, which describes the distance between the highest point of projection and the caudal implant edge, establishes a reference parameter for the correct placement of the areola in relation to the inframammary fold.

Information on Material

POLYTECH Health & Aesthetics implants are filled with the highly cross-linked EasyFit™ gel, which stands out due to its natural consistency. Because of the gel's "memory effect", it always returns to its original shape after moderate mechanical deformation.

The chemically and mechanically resistant, elastic shell of the implants consists of several layers of a silicone elastomer. In order to maintain a high level of safety, the shell of the silicone gel-filled implants is equipped with an additional diffusion barrier layer.

The surface of the implant shell is Microthane®-covered (micropolyurethane foam). This surface is known for its reliability and has been proven to contribute to a reduction of capsular contracture rates¹.

Production of the implants is based on a quality management according to ISO 9001 and ISO 13485. All materials used are of implantable medical quality.

1: (Handel, 1991; Gasperoni, 1992; Pennisi, 1990; Shapiro, 1989; Baudelot, 1989; Artz, 1988; Hermann, 1984; Eyssen, 1984; Schatten, 1984; Hester et al., 2001; Vázquez, 2007) (Handel, 2006)

Indications

- Breast reconstruction after mastectomy or ablatio mammae
- Contour augmentation and correction
- Adjustment of asymmetries, aplasia or hypoplasia
- Combined breast and thoracic deformities
- Ptosis or atrophy
- Implant exchange

Important Information

The implants are delivered individually packed and sterile.

Attention: Please note the detailed instructions for use included with each product.



POLYTECH Health & Aesthetics GmbH

Altheimer Str. 32 • 64807 Dieburg • Germany

+49 (0) 60 71 98 63 0 • +49 (0) 60 71 98 63 30

info@polytechhealth.com • www.polytech-health-aesthetics.com