



Mème® LP/S



Mème® MP/S



Mème® HP/S



Mème® XP/S

Mème® *SublimeLine*®

POLYsmooth™

Silicone-gel-filled mammary implants, with smooth surface, round base, central profile, maximum shape control

Implants made by POLYTECH – QUALITY made in Germany



POLYTECH
health & aesthetics 

Même® POLYsmooth™																	
Base Width (mm)		LP Low Projection			MP Moderate Projection			HP High Projection			XP Extra-High Projection						
A	C	Order n°	Projection (mm)		Order n°	Projection (mm)		Order n°	Projection (mm)		Order n°	Projection (mm)					
			B	D		B	D		B	D		B	D				
				Vol. (ml)			Vol. (ml)			Vol. (ml)			Vol. (ml)				
70	70	10724-055*	23	43	55	10725-070*	28	47	70	10726-090*	34	51	90	10727-105*	40	55	105
75	75	10724-065 ¹	24	46	65	10725-090 ¹	30	50	90	10726-110 ¹	36	54	110	10727-125	43	59	125
80	80	10724-085*	24	49	85	10725-105*	31	53	105	10726-125*	37	58	125	10727-155*	45	62	155
85	85	10724-095 ¹	25	51	95	10725-120	32	56	120	10726-155	39	61	155	10727-180	47	66	180
90	90	10724-110*	26	54	110	10725-145*	33	59	145	10726-180*	40	64	180	10727-215*	49	69	215
95	95	10724-125	27	57	125	10725-165	35	62	165	10726-210	42	67	210	10727-250	51	73	250
100	100	10724-150*	28	60	150	10725-195*	36	65	195	10726-235*	44	71	235	10727-285*	53	77	285
105	105	10724-165	29	62	165	10725-220	37	68	220	10726-275	45	74	275	10727-330	55	80	330
110	110	10724-185*	29	65	185	10725-255*	38	71	255	10726-315*	47	77	315	10727-380*	57	84	380
115	115	10724-215	30	68	215	10725-285	40	74	285	10726-360	49	80	360	10727-430	59	87	430
120	120	10724-235*	31	71	235	10725-320*	41	77	320	10726-400*	51	83	400	10727-485*	61	91	485
125	125	10724-270	32	73	270	10725-360	43	80	360	10726-450	53	87	450	10727-545	63	94	545
130	130	10724-295*	33	76	295	10725-400*	44	83	400	10726-505*	55	90	505	10727-610*	65	98	610
135	135	10724-330	34	79	330	10725-445	45	86	445	10726-560	57	93	560	10727-675	67	101	675
140	140	10724-365*	35	82	365	10725-495*	47	89	495	10726-620*	59	96	620	10727-755*	69	105	755
145	145	10724-400 ¹	35	85	400	10725-545 ¹	48	92	545	10726-690 ¹	61	100	690	10727-830 ¹	71	108	830
150	150	10724-440*	36	87	440	10725-600*	49	95	600	10726-760*	62	103	760	10727-920*	73	112	920
155	155	10724-485 ¹	37	90	485	10725-660 ¹	51	98	660	10726-835 ¹	64	106	835	10727-1010 ¹	75	115	1010

Approximate measurements. Sizes marked 1 and set in grey as well as further sizes will be produced on demand. Silicone gel-filled, smooth sizers (non sterile, re-sterilisable) are available for products marked with an asterisk *.

Product Description

The silicone-gel-filled Même® *SublimeLine*® implants with POLYsmooth™ surface have a round base and a central projection. The filling level of the *SublimeLine*® implants is maximized to provide more volume.

The consistency of these high-quality mammary implants is very close to the haptics and motional properties of the mammae. The implants are available in four projections: Low Profile (LP), Moderate Profile (MP), High Profile (HP), and Extra High Profile (XP).

The special curve measurement D, which describes the distance between the highest point of projection and the caudal implant edge, establishes a reference parameter for the correct placement of the areola in relation to the inframammary fold.

Information on Material

POLYTECH Health & Aesthetics implants are filled with the highly cross-linked EasyFit™ gel, which stands out due to its natural consistency. Because of the gel's "memory effect", it always returns to its original shape after moderate mechanical deformation.

The chemically and mechanically resistant, elastic shell of the implants consists of several layers of a silicone elastomer. In order to maintain a high level of safety, the shell of the silicone gel-filled implants is equipped with an additional diffusion barrier layer. The implant shell has a POLYsmooth™ smooth surface.

Production of the implants is based on a quality management according to ISO 9001 and ISO 13485. All materials used are of implantable medical quality.

Indications

- Breast reconstruction after mastectomy or ablatio mammae
- Contour augmentation and correction
- Adjustment of asymmetries, aplasia or hypoplasia
- Combined breast and thoracic deformities
- Ptosis or atrophy
- Implant exchange

Important Information

The implants are delivered individually packed and sterile.

Attention: Please note the detailed instructions for use included with each product.



POLYTECH Health & Aesthetics GmbH

Altheimer Str. 32 • 64807 Dieburg • Germany

+49 (0) 60 71 98 63 0 • +49 (0) 60 71 98 63 30

info@polytechhealth.com • www.polytech-health-aesthetics.com