



McKinley Medical L.L.P.
 4080 Youngfield Street
 Wheat Ridge, CO 80033 U.S.A
 (303) 420.9589
 (303) 420.8585 Fax

ACCUFUSER™/Plus Pain Pump

Quick Reference - Daily Delivery Chart (mg) of Local Anesthetics

Drug Used	Max. Daily mg	Concentration		Infusion Rate (ml/hr)		
		%	mg/ml	2	4	5
Bupivacaine	400	0.25%	2.5	120	240	300
		0.50%	5	240	(480)	(600)
Lidocaine	300	0.50%	5	240	(480)	(600)
Lidocaine with Epinephrine	500	0.50%	5	240	480	(600)
		1.00%	10	480	(960)	(1200)
Levobupivacaine (Chirocaine)	695	0.25%	2.5	120	240	300
		0.50%	5	240	480	600
Ropivacaine	770	0.20%	2	96	192	240
		0.50%	5	240	480	600
		0.75%	7.5	360	720	(900)
		1.00%	10	480	(960)	(1200)

Values in () exceed maximum recommended daily dose.