


**Effective Pain Control
Using Less Medication**

Clinical research shows that pumps providing continuous basal infusion as well as patient-controlled bolus dosing provide equivalent or superior analgesia with less medication consumption compared to continuous infusion or bolus-dosing delivered alone.* ACCUFUSER is a non-electronic balloon infuser that provides you the therapeutic advantage of combined continuous infusion plus patient-controlled bolus-dosing.



**Select Your
Therapeutic Options**

5 Basal Infusion Rates with Bolus:

- 1.0 ml/hr 5.0 ml/hr
- 2.0 ml/hr 8.0 ml/hr
- 4.0 ml/hr

3 Bolus Dose Volumes:

- 0.5 ml 2.0 ml
- 1.0 ml

3 Bolus Dose Lockout Intervals:

- 15 min. 60 min.
- 30 min.

4 Medication Reservoir Sizes:

- 60 ml 275 ml
- 100 ml 550 ml

Select the ACCUFUSER That's Right for You

ACCUFUSER Model	Reservoir Volume	Basal Rate	Bolus Volume	Lockout Interval	Max Bolus Per Hour	Max Hourly Volume
M1015M	100 ml	1 ml/hr	1 ml	15 min.	4 ml	5 ml
M2015M	100 ml	2 ml/hr	1 ml	15 min.	4 ml	6 ml
T2015M	100 ml	2 ml/hr	2 ml	15 min.	8 ml	10 ml
M2060M	100 ml	2 ml/hr	1 ml	60 min.	1 ml	3 ml
T2060M	100 ml	2 ml/hr	2 ml	60 min.	2 ml	4 ml
T2015L	275 ml	2 ml/hr	2 ml	15 min.	8 ml	10 ml
T2060L	275 ml	2 ml/hr	2 ml	60 min.	2 ml	4 ml
T4030L	275 ml	4 ml/hr	2 ml	30 min.	4 ml	8 ml
T4060L	275 ml	4 ml/hr	2 ml	60 min.	2 ml	6 ml
T5015L	275 ml	5 ml/hr	2 ml	15 min.	8 ml	13 ml
T5060L	275 ml	5 ml/hr	2 ml	60 min.	2 ml	7 ml
T8015L	275 ml	8 ml/hr	2 ml	15 min.	8 ml	16 ml
T8060L	275 ml	8 ml/hr	2 ml	60 min.	2 ml	10 ml
T5015XL	550 ml	5 ml/hr	2 ml	15 min.	8 ml	13 ml
T5060XL	550 ml	5 ml/hr	2 ml	60 min.	2 ml	7 ml
T8015XL	550 ml	8 ml/hr	2 ml	15 min.	8 ml	16 ml
T8060XL	550 ml	8 ml/hr	2 ml	60 min.	2 ml	10 ml

Note: These are the most commonly used ACCUFUSER models. Ask your McKinley Medical representative about additional models to meet your needs.

* Singelyn FJ, Seguy S, Gouverneur JM. Interscalene brachial plexus analgesia after open shoulder surgery: continuous versus patient-controlled infusion. *Anesth Analg* 1999; 89 (5): 1216 – 20.



McKinley Medical

McKinley Medical, LLC
4080 Youngfield Street
Wheat Ridge, CO 80033 USA
303.420.9569
303.420.8585 (fax)
www.controlpostoppain.com

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

©2005 by McKinley Medical, LLC. All rights reserved. ACCUFUSER and the McKinley Medical mountain logo are registered trademarks, and The Beauty of Personal Post-Op Pain Control and McKinley Saturation Catheter are trademarks of McKinley Medical, LLC. U.S. Patents: 6,024,724; 6,689,110B2; 6,676,643B2. **PN 203656 Rev. AA 01/2005**