



Mème® LP



Mème® MP



Mème® HP



Mème® XP

Mème® *SublimeLine®*

MESMO® sensitive

Silicone-gel-filled mammary implants, with micro-textured surface, round base, central profile, maximum shape control

Implants made by POLYTECH – QUALITY made in Germany



POLYTECH
health & aesthetics 

Même® – micro-textured													
		MP Moderate Projection			HP High Projection			XP Extra-High Projection					
Base (mm)		Order n°	Projection (mm)		Vol. (ml)	Order n°	Projection (mm)		Vol. (ml)	Order n°	Projection (mm)		Vol. (ml)
A	B		C	D		C	D				C	D	
70	70	15725-070 ¹ *	28	47	70	15726-090 ¹ *	34	51	90				
75	75	15725-090 ¹	30	50	90	15726-110 ¹	36	54	110				
80	80	15725-105 ¹ *	31	56	105	15726-125 ¹ *	37	58	125				
85	85	15725-120 ¹	32	56	120	15726-155 ¹	39	61	155				
90	90	15725-145*	33	59	145	15726-180*	40	64	180	15727-215 ¹ *	49	69	215
95	95	15725-165	35	62	165	15726-210	42	67	210	15727-250 ¹	51	73	250
100	100	15725-195*	36	65	195	15726-235*	44	71	235	15727-285 ¹ *	53	77	285
103	103					15726-255	45	72	255	15727-305 ¹	54	78	305
105	105	15725-220	37	68	220	15726-275	45	74	275	15727-330	55	80	330
108	108					15726-295	46	75	295	15727-355*	56	82	355
110	110	15725-255*	38	71	255	15726-315*	47	77	315	15727-380*	57	84	380
113	113					15726-335	48	79	335				
115	115	15725-285	40	74	285	15726-360	49	80	360	15727-430	59	87	430
118	118	15725-300	40	75	300								
120	120	15725-320*	41	77	320	15726-400*	51	83	400	15727-485 ¹ *	61	91	485
125	125	15725-360	43	80	360	15726-450	53	87	450				
130	130	15725-400*	44	83	400	15726-505*	55	90	505				
135	135	15725-445	45	86	445	15726-560	57	93	560				
140	140	15725-495 ¹ *	47	89	495	15726-620 ¹ *	59	96	620				
145	145	15725-545 ¹	48	92	545	15726-690 ¹	61	100	690				
150	150	15725-600 ¹ *	49	95	600	15726-760 ¹ *	62	103	760				
155	155	15725-660 ¹	51	98	660	15726-835 ¹	64	106	835				

Approximate measurements. Sizes highlighted in violet and marked (1) will be produced on demand; supply = 6 weeks after order. Silicone gel-filled, smooth sizers (sterile, re-sterilisable) are available for products marked with an asterisk *. Further sizes will be produced according to our terms for custom-made products.

Product Description

The silicone-gel-filled Même® *SublimeLine*® implants with MESMO®sensitive surface have a round base and a central projection. The filling level of the *SublimeLine*® implants is maximized to provide more volume.

The consistency of these high-quality mammary implants is very close to the haptics and motional properties of the mammae. The implants are available in four projections: Low Profile (LP), Moderate Profile (MP), High Profile (HP), and Extra High Profile (XP).

The special curve measurement D, which describes the distance between the highest point of projection and the caudal implant edge, establishes a reference parameter for the correct placement of the areola in relation to the inframammary fold.

Information on Material

POLYTECH Health & Aesthetics implants are filled with the highly cross-linked EasyFit™ gel, which stands out due to its natural consistency. Because of the gel's "memory effect", it always returns to its original shape after moderate mechanical deformation.

The chemically and mechanically resistant, elastic shell of the implants consists of several layers of a silicone elastomer. In order to maintain a high level of safety, the shell of the silicone gel-filled implants is equipped with an additional diffusion barrier layer. The surface of the implant shell is MESMO®sensitive micro-textured.

Production of the implants is based on a quality management according to ISO 9001 and ISO 13485. All materials used are of implantable medical quality.

Indications

- Breast reconstruction after mastectomy or ablatio mammae
- Contour augmentation and correction
- Adjustment of asymmetries, aplasia or hypoplasia
- Combined breast and thoracic deformities
- Ptosis or atrophy
- Implant exchange

Important Information

The implants are delivered individually packed and sterile.

Attention: Please note the detailed instructions for use included with each product.

POLYTECH Health & Aesthetics GmbH

Altheimer Str. 32 • 64807 Dieburg • Germany

+49 (0) 60 71 98 63 0 • +49 (0) 60 71 98 63 30

info@polytechhealth.com • www.polytech-health-aesthetics.com

