



Mème® LP/T



Mème® MP/T



Mème® HP/T



Mème® XP/T

Mème® *SublimeLine*®

POLYtxt®

Silicone-gel-filled mammary implants, with textured surface, round base, central profile, maximum shape control

Implants made by POLYTECH – QUALITY made in Germany



Même® POLYtxt®																	
Base Width (mm)		LP Low Projection			MP Moderate Projection			HP High Projection			XP Extra-High Projection						
A	C	Order n°	Projection (mm)		Order n°	Projection (mm)		Order n°	Projection (mm)		Order n°	Projection (mm)					
			B	D		B	D		B	D		B	D				
				Vol. (ml)			Vol. (ml)			Vol. (ml)			Vol. (ml)				
70	70	20724-055*	23	43	55	20725-070*	28	47	70	20726-090*	34	51	90	20727-105*	40	55	105
75	75	20724-065 ¹	24	46	65	20725-090 ¹	30	50	90	20726-110 ¹	36	54	110	20727-125	43	59	125
80	80	20724-085*	24	49	85	20725-105*	31	53	105	20726-125*	37	58	125	20727-155*	45	62	155
85	85	20724-095 ¹	25	51	95	20725-120	32	56	120	20726-155	39	61	155	20727-180	47	66	180
90	90	20724-110*	26	54	110	20725-145*	33	59	145	20726-180*	40	64	180	20727-215*	49	69	215
95	95	20724-125	27	57	125	20725-165	35	62	165	20726-210	42	67	210	20727-250	51	73	250
100	100	20724-150*	28	60	150	20725-195*	36	65	195	20726-235*	44	71	235	20727-285*	53	77	285
105	105	20724-165	29	62	165	20725-220	37	68	220	20726-275	45	74	275	20727-330	55	80	330
110	110	20724-185*	29	65	185	20725-255*	38	71	255	20726-315*	47	77	315	20727-380*	57	84	380
115	115	20724-215	30	68	215	20725-285	40	74	285	20726-360	49	80	360	20727-430	59	87	430
120	120	20724-235*	31	71	235	20725-320*	41	77	320	20726-400*	51	83	400	20727-485*	61	91	485
125	125	20724-270	32	73	270	20725-360	43	80	360	20726-450	53	87	450	20727-545	63	94	545
130	130	20724-295*	33	76	295	20725-400*	44	83	400	20726-505*	55	90	505	20727-610*	65	98	610
135	135	20724-330	34	79	330	20725-445	45	86	445	20726-560	57	93	560	20727-675	67	101	675
140	140	20724-365*	35	82	365	20725-495*	47	89	495	20726-620*	59	96	620	20727-755*	69	105	755
145	145	20724-400 ¹	35	85	400	20725-545 ¹	48	92	545	20726-690 ¹	61	100	690	20727-830 ¹	71	108	830
150	150	20724-440*	36	87	440	20725-600*	49	95	600	20726-760*	62	103	760	20727-920*	73	112	920
155	155	20724-485 ¹	37	90	485	20725-660 ¹	51	98	660	20726-835 ¹	64	106	835	20727-1010 ¹	75	115	1010

Approximate measurements. Sizes marked 1 and set in grey as well as further sizes will be produced on demand. Silicone gel-filled, smooth sizers (non sterile, re-sterilisable) are available for products marked with an asterisk *.

Product Description

The silicone-gel-filled Même® *SublimeLine*® implants with POLYtxt® surface have a round base and a central projection. The filling level of the *SublimeLine*® implants is maximized to provide more volume.

The consistency of these high-quality mammary implants is very close to the haptics and motional properties of the mammae. The implants are available in four projections: Low Profile (LP), Moderate Profile (MP), High Profile (HP), and Extra High Profile (XP).

The special curve measurement D, which describes the distance between the highest point of projection and the caudal implant edge, establishes a reference parameter for the correct placement of the areola in relation to the inframammary fold.

Information on Material

POLYTECH Health & Aesthetics implants are filled with the highly cross-linked EasyFit™ gel, which stands out due to its natural consistency. Because of the gel's "memory effect", it always returns to its original shape after moderate mechanical deformation.

The chemically and mechanically resistant, elastic shell of the implants consists of several layers of a silicone elastomer. In order to maintain a high level of safety, the shell of the silicone gel-filled implants is equipped with an additional diffusion barrier layer.

The surface of the implant shell is POLYtxt® textured. This surface is known for its reliability and has been proven to contribute to a reduction of capsular contracture rates¹.

Production of the implants is based on a quality management according to ISO 9001 and ISO 13485. All materials used are of implantable medical quality.

¹: textured implants contracture rates: 8–15% (Iwuagwu + Frame, 1997), which is 15% less than with smooth implants (Handel, 2006)

Indications

- Breast reconstruction after mastectomy or ablatio mammae
- Contour augmentation and correction
- Adjustment of asymmetries, aplasia or hypoplasia
- Combined breast and thoracic deformities
- Ptosis or atrophy
- Implant exchange

Important Information

The implants are delivered individually packed and sterile.

Attention: Please note the detailed instructions for use included with each product.



POLYTECH Health & Aesthetics GmbH

Altheimer Str. 32 • 64807 Dieburg • Germany

☎ +49 (0) 60 71 98 63 0 • 📠 +49 (0) 60 71 98 63 30

✉ info@polytechhealth.com • www.polytech-health-aesthetics.com