



Opticon® LP/M



Opticon® MP/M



Opticon® HP/M



Opticon® XP/M

Opticon® *SublimeLine*®

Microthane®

Silicone-gel-filled mammary implants, with Microthane® surface, short oval base, anatomical profile, maximum shape control

Implants made by POLYTECH – QUALITY made in Germany



Opticon® Microthane®																	
Base Width (mm)		LP Low Projection			MP Moderate Projection			HP High Projection			XP Extra-High Projection						
A	C	Order n°	Projection (mm)		Vol. (ml)	Order n°	Projection (mm)		Vol. (ml)	Order n°	Projection (mm)		Vol. (ml)				
73	62	30744-045*	25	32	45	30745-055*	30	39	55	30746-065*	36	44	65	30747-075*	42	48	75
78	66	30744-050 ¹	26	34	50	30745-065 ¹	32	41	65	30746-075 ¹	38	46	75	30747-095 ¹	45	51	95
83	70	30744-060*	26	35	60	30745-075*	33	43	75	30746-095*	39	48	95	30747-115*	47	54	115
88	74	30744-070 ¹	27	37	70	30745-095 ¹	34	45	95	30746-115	41	51	115	30747-130	49	56	130
93	79	30744-085*	28	39	85	30745-105*	35	47	105	30746-130*	42	53	130	30747-160*	51	59	160
98	83	30744-095 ¹	29	41	95	30745-125	37	49	125	30746-155	44	55	155	30747-185	53	62	185
103	87	30744-110*	30	42	110	30745-145*	38	51	145	30746-175*	46	58	175	30747-215*	55	65	215
108	91	30744-120	31	44	120	30745-165	39	53	165	30746-205	47	60	205	30747-240	57	67	240
113	95	30744-140*	31	46	140	30745-185*	40	55	185	30746-230*	49	62	230	30747-280*	59	70	280
118	100	30744-160	32	48	160	30745-210	42	57	210	30746-265	51	65	265	30747-320	61	73	320
123	104	30744-175*	33	49	175	30745-235*	43	59	235	30746-295*	53	67	295	30747-360*	63	75	360
128	108	30744-200	34	51	200	30745-265	45	61	265	30746-335	55	69	335	30747-400	65	78	400
133	112	30744-220*	35	53	220	30745-295*	46	63	295	30746-375*	57	72	375	30747-450*	67	81	450
138	116	30744-240	36	55	240	30745-330	47	65	330	30746-420	59	74	420	30747-500	69	83	500
143	121	30744-270*	37	56	270	30745-365*	49	67	365	30746-460*	61	76	460	30747-560*	71	86	560
148	125	30744-295	37	58	295	30745-400	50	69	400	30746-510	63	79	510	30747-615	73	89	615
153	129	30744-325*	38	60	325	30745-445*	51	71	445	30746-560*	64	81	560	30747-680*	75	91	680
158	133	30744-360	39	62	360	30745-490	53	73	490	30746-615	66	84	615	30747-750 ¹	77	94	750

Approximate measurements. Sizes marked 1 and set in grey as well as further sizes will be produced on demand. Silicone gel-filled, smooth sizers (non sterile, re-sterilisable) are available for products marked with an asterisk *.

Product Description

The silicone-gel-filled Opticon® *SublimeLine*® implants with Microthane® surface have a short oval base and an anatomical projection with the highest point in the lower third. The filling level of the *SublimeLine*® implants is maximized to provide more volume.

The consistency of these high-quality mammary implants is very close to the haptics and motional properties of the mammae. The implants are available in four projections: Low Profile (LP), Moderate Profile (MP), High Profile (HP), and Extra High Profile (XP).

The special curve measurement D, which describes the distance between the highest point of projection and the caudal implant edge, establishes a reference parameter for the correct placement of the areola in relation to the inframammary fold.

Information on Material

POLYTECH Health & Aesthetics implants are filled with the highly cross-linked EasyFit™ gel, which stands out due to its natural consistency. Because of the gel's "memory effect", it always returns to its original shape after moderate mechanical deformation.

The chemically and mechanically resistant, elastic shell of the implants consists of several layers of a silicone elastomer. In order to maintain a high level of safety, the shell of the silicone gel-filled implants is equipped with an additional diffusion barrier layer.

The surface of the implant shell is Microthane®-covered (micropolyurethane foam). This surface is known for its reliability and has been proven to contribute to a reduction of capsular contracture rates¹.

Production of the implants is based on a quality management according to ISO 9001 and ISO 13485. All materials used are of implantable medical quality.

1: (Handel, 1991; Gasperoni, 1992; Pennisi, 1990; Shapiro, 1989; Baudelot, 1989; Artz, 1988; Hermann, 1984; Eyssen, 1984; Schatten, 1984; Hester et al., 2001; Vázquez, 2007) (Handel, 2006)

Indications

- Breast reconstruction after mastectomy or ablatio mammae
- Contour augmentation and correction
- Adjustment of asymmetries, aplasia or hypoplasia
- Combined breast and thoracic deformities
- Ptosis or atrophy
- Implant exchange

Important Information

The implants are delivered individually packed and sterile.

Attention: Please note the detailed instructions for use included with each product.



POLYTECH Health & Aesthetics GmbH

Altheimer Str. 32 • 64807 Dieburg • Germany

+49 (0) 60 71 98 63 0 • +49 (0) 60 71 98 63 30

info@polytechhealth.com • www.polytech-health-aesthetics.com