



Opticon® LP/M



Opticon® MP/M



Opticon® HP/M



Opticon® XP/M

Opticon® *SublimeLine*®

POLYtxt®

Silicone-gel-filled mammary implants, with textured surface, short oval base, anatomical profile, maximum shape control

Implants made by POLYTECH – QUALITY made in Germany



POLYTECH
health & aesthetics 

Opticon® – POLYtxt®													
Base (mm)		LP Low Projection			MP Moderate Projection			HP High Projection			XP Extra-High Projection		
A	B	Order n°	Projection (mm)	Vol. (ml)	Order n°	Projection (mm)	Vol. (ml)	Order n°	Projection (mm)	Vol. (ml)	Order n°	Projection (mm)	Vol. (ml)
70	59	20744-045 ^{1*}	23 30	45	20745-055 ^{1*}	28 35	55	20746-065 ^{1*}	34 40	65	20747-075 ^{1*}	40 44	75
75	63	20744-050 ¹	24 32	50	20745-065 ¹	30 37	65	20746-075 ¹	36 42	75	20747-095 ¹	43 47	95
80	67	20744-060 ^{1*}	24 34	60	20745-075 ^{1*}	31 39	75	20746-095 ^{1*}	37 44	95	20747-115 ^{1*}	45 50	115
85	71	20744-070 ¹	25 35	70	20745-095 ¹	32 41	95	20746-115	39 47	115	20747-130	47 52	130
90	76	20744-085 ^{1*}	26 37	85	20745-105 ^{1*}	33 43	105	20746-130 ^{1*}	40 49	130	20747-160 ^{1*}	49 55	160
95	80	20744-095 ¹	27 39	95	20745-125	35 45	125	20746-155	42 51	155	20747-185	51 58	185
100	84	20744-110 ^{1*}	28 41	110	20745-145 ^{1*}	36 47	145	20746-175 ^{1*}	44 54	175	20747-215 ^{1*}	53 61	215
105	88	20744-120	29 42	120	20745-165	37 49	165	20746-205	45 56	205	20747-240	55 63	240
110	92	20744-140 ^{1*}	29 44	140	20745-185 ^{1*}	38 51	185	20746-230 ^{1*}	47 58	230	20747-280 ^{1*}	57 66	280
115	97	20744-160	30 46	160	20745-210	40 53	210	20746-265	49 61	265	20747-320	59 69	320
120	101	20744-175 ^{1*}	31 48	175	20745-235 ^{1*}	41 55	235	20746-295 ^{1*}	51 63	295	20747-360 ^{1*}	61 71	360
125	105	20744-200	32 49	200	20745-265	43 57	265	20746-335	53 65	335	20747-400	63 74	400
130	109	20744-220 ^{1*}	33 51	220	20745-295 ^{1*}	44 59	295	20746-375 ^{1*}	55 68	375	20747-450 ^{1*}	65 77	450
135	113	20744-240	34 53	240	20745-330	45 61	330	20746-420	57 70	420	20747-500	67 79	500
140	118	20744-270 ^{1*}	35 55	270	20745-365 ^{1*}	47 63	365	20746-460 ^{1*}	59 72	460	20747-560 ^{1*}	69 82	560
145	122	20744-295	35 56	295	20745-400	48 65	400	20746-510	61 75	510	20747-615	71 85	615
150	126	20744-325 ^{1*}	36 58	325	20745-445 ^{1*}	49 67	445	20746-560 ^{1*}	62 77	560	20747-680 ^{1*}	73 87	680
155	130	20744-360	37 60	360	20745-490	51 69	490	20746-615	64 80	615	20747-750 ¹	75 90	750

Approximate measurements. Sizes highlighted in color and marked (1) will be produced on demand; supply ≈ 6 weeks after order. Silicone gel-filled, smooth sizers (sterile, re-sterilisable) are available for products marked with an asterisk (*). Further sizes will be produced according to our terms for custom-made products.

Product Description

The silicone-gel-filled Opticon® *SublimeLine*® implants with POLYtxt® surface have a short oval base and an anatomical projection with the highest point in the lower third. The filling level of the *SublimeLine*® implants is maximized to provide more volume.

The consistency of these high-quality mammary implants is very close to the haptics and motional properties of the mammae. The implants are available in four projections: Low Profile (LP), Moderate Profile (MP), High Profile (HP), and Extra High Profile (XP).

The special curve measurement D, which describes the distance between the highest point of projection and the caudal implant edge, establishes a reference parameter for the correct placement of the areola in relation to the inframammary fold.

Information on Material

POLYTECH Health & Aesthetics implants are filled with the highly cross-linked EasyFit™ gel, which stands out due to its natural consistency. Because of the gel's "memory effect", it always returns to its original shape after moderate mechanical deformation.

The chemically and mechanically resistant, elastic shell of the implants consists of several layers of a silicone elastomer. In order to maintain a high level of safety, the

shell of the silicone gel-filled implants is equipped with an additional diffusion barrier layer.

The surface of the implant shell is POLYtxt® textured. This surface is known for its reliability and has been proven to contribute to a reduction of capsular contracture rates¹.

Production of the implants is based on a quality management according to ISO 9001 and ISO 13485. All materials used are of implantable medical quality.

¹: textured implants contracture rates: 8–15% (Iwuagwu + Frame, 1997), which is 15% less than with smooth implants (Handel, 2006)

Indications

- Breast reconstruction after mastectomy or ablatio mammae
- Contour augmentation and correction
- Adjustment of asymmetries, aplasia or hypoplasia
- Combined breast and thoracic deformities
- Ptosis or atrophy
- Implant exchange

Important Information

The implants are delivered individually packed and sterile.

Attention: Please note the detailed instructions for use included with each product.



POLYTECH Health & Aesthetics GmbH

Altheimer Str. 32 • 64807 Dieburg • Germany

+49 (0) 60 71 98 63 0 • +49 (0) 60 71 98 63 30

info@polytechhealth.com • www.polytech-health-aesthetics.com