



Replicon® LP



Replicon® MP



Replicon® HP



Replicon® XP

Replicon® *SublimeLine*®

MESMO® sensitive

Silicone-gel-filled mammary implants,
micro-textured surface, round base, anatomical
profile, maximum shape control

Implants made by POLYTECH – QUALITY made in Germany



POLYTECH
health & aesthetics 

Replicon® – micro-textured													
Base (mm)		MP Moderate Projection			HP High Projection			XP Extra-High Projection					
A	B	Order n°	Projection (mm)		Vol. (ml)	Order n°	Projection (mm)		Vol. (ml)	Order n°	Projection (mm)		Vol. (ml)
			C	D			C	D			C	D	
100	100					15736-210*	44	57	210	15737-255 ¹ *	53	63	255
105	105					15736-240	45	59	240	15737-290	55	66	290
110	110	15735-220*	38	55	220	15736-275*	47	62	275	15737-330*	57	69	330
115	115	15735-255	40	57	255	15736-315	49	64	315	15737-380	59	72	380
120	120	15735-280*	41	59	280	15736-350*	51	67	350	15737-425*	61	75	425
125	125	15735-320	43	61	320	15736-395	53	69	395	15737-480	63	78	480
130	130	15735-350*	44	63	350	15736-445*	55	72	445	15737-535 ¹ *	65	80	535
135	135	15735-390	45	66	390	15736-495	57	74	495	15737-600 ¹	67	83	600
140	140	15735-435*	47	68	435	15736-550*	59	77	550				
145	145	15735-480 ¹	48	70	480	15736-605 ¹	61	79	605				
150	150	15735-530 ¹ *	49	72	530	15736-670 ¹ *	62	82	670				
155	155					15736-735 ¹ *	64	84	735				

Approximate measurements. Sizes highlighted in violet and marked (1) will be produced on demand; supply ≈ 6 weeks after order. Silicone gel-filled, smooth sizers (sterile, re-sterilisable) are available for products marked with an asterisk *. Further sizes will be produced according to our terms for custom-made products.

Product Description

The silicone-gel-filled Replicon® *SublimeLine*® implants with MESMO® sensitive surface have a round base and an anatomical projection with the highest point in the lower half. The filling level of the *SublimeLine*® implants is maximized to provide more volume.

The consistency of these high-quality mammary implants is very close to the haptics and motional properties of the mammae. The implants are available in four projections: Low Profile (LP), Moderate Profile (MP), High Profile (HP), and Extra High Profile (XP).

The special curve measurement D, which describes the distance between the highest point of projection and the caudal implant edge, establishes a reference parameter for the correct placement of the areola in relation to the inframammary fold.

Information on Material

POLYTECH Health & Aesthetics implants are filled with the highly cross-linked EasyFit™ gel, which stands out due to its natural consistency. Because of the gel's "memory effect", it always returns to its original shape after moderate mechanical deformation.

The chemically and mechanically resistant, elastic shell of the implants consists of several layers of a silicone elastomer. In order to maintain a high level of safety, the

shell of the silicone gel-filled implants is equipped with an additional diffusion barrier layer. The surface of the implant shell is MESMO® sensitive micro-textured.

Production of the implants is based on a quality management according to ISO 9001 and ISO 13485. All materials used are of implantable medical quality.

Indications

- Breast reconstruction after mastectomy or ablatio mammae
- Contour augmentation and correction
- Adjustment of asymmetries, aplasia or hypoplasia
- Combined breast and thoracic deformities
- Ptosis or atrophy
- Implant exchange

Important Information

The implants are delivered individually packed and sterile.

Attention: Please note the detailed instructions for use included with each product.



POLYTECH Health & Aesthetics GmbH

Altheimer Str. 32 • 64807 Dieburg • Germany

+49 (0) 60 71 98 63 0 • +49 (0) 60 71 98 63 30

info@polytechhealth.com • www.polytech-health-aesthetics.com