



Replicon® LP/T



Replicon® MP/T



Replicon® HP/T



Replicon® XP/T

Replicon® *SublimeLine*®

POLYtxt®

Silicone-gel-filled mammary implants,
with textured surface, round base, anatomical
profile, maximum shape control

Implants made by POLYTECH – QUALITY made in Germany



POLYTECH
health & aesthetics 

Replicon® POLYtxt®																	
Base with (mm)		LP Low Projection			MP Moderate Projection			HP High Projection			XP Extra-High Projection						
A	C	Order n°	Projection (mm)		Vol. (ml)	Order n°	Projection (mm)		Vol. (ml)	Order n°	Projection (mm)		Vol. (ml)				
			B	D			B	D			B	D					
70	70	20734-050*	23	33	50	20735-065*	28	37	65	20736-075*	34	42	75	20737-095*	40	46	95
75	75	20734-060 ¹	24	34	60	20735-075 ¹	30	39	75	20736-095 ¹	36	44	95	20737-110	43	49	110
80	80	20734-070*	24	36	70	20735-095*	31	41	95	20736-115*	37	47	115	20737-130*	45	52	130
85	85	20734-085 ¹	25	38	85	20735-110	32	44	110	20736-130	39	49	130	20737-160	47	55	160
90	90	20734-100*	26	40	100	20735-125*	33	46	125	20736-160*	40	52	160	20737-185*	49	58	185
95	95	20734-110	27	42	110	20735-150	35	48	150	20736-180	42	54	180	20737-220	51	61	220
100	100	20734-125*	28	44	125	20735-170*	36	50	170	20736-210*	44	57	210	20737-255*	53	63	255
105	105	20734-150	29	46	150	20735-195	37	52	195	20736-240	45	59	240	20737-290	55	66	290
110	110	20734-165*	29	48	165	20735-220*	38	55	220	20736-275*	47	62	275	20737-330*	57	69	330
115	115	20734-185	30	50	185	20735-255	40	57	255	20736-315	49	64	315	20737-380	59	72	380
120	120	20734-210*	31	52	210	20735-280*	41	59	280	20736-350*	51	67	350	20737-425*	61	75	425
125	125	20734-235	32	54	235	20735-320	43	61	320	20736-395	53	69	395	20737-480	63	78	480
130	130	20734-265*	33	56	265	20735-350*	44	63	350	20736-445*	55	72	445	20737-535*	65	80	535
135	135	20734-290	34	58	290	20735-390	45	66	390	20736-495	57	74	495	20737-600	67	83	600
140	140	20734-350	35	59	320	20735-435*	47	68	435	20736-550*	59	77	550	20737-665*	69	86	665
145	145	20734-320	35	61	350	20735-480	48	70	480	20736-605 ¹	61	79	605	20737-730 ¹	71	89	730
150	150	20734-390*	36	63	390	20735-530*	49	72	530	20736-670*	62	82	670	20737-810*	73	92	810
155	155	20734-425	37	65	425	20735-580	51	74	580	20736-735 ¹	64	84	735	20737-890 ¹	75	95	890

Approximate measurements. Sizes marked 1 and set in grey as well as further sizes will be produced on demand. Silicone gel-filled, smooth sizers (non sterile, re-sterilisable) are available for products marked with an asterisk *.

Product Description

The silicone-gel-filled Replicon® *SublimeLine*® implants with POLYtxt® surface have a round base and an anatomical projection with the highest point in the lower half. The filling level of the *SublimeLine*® implants is maximized to provide more volume.

The consistency of these high-quality mammary implants is very close to the haptics and motional properties of the mammae. The implants are available in four projections: Low Profile (LP), Moderate Profile (MP), High Profile (HP), and Extra High Profile (XP).

The special curve measurement D, which describes the distance between the highest point of projection and the caudal implant edge, establishes a reference parameter for the correct placement of the areola in relation to the inframammary fold.

Information on Material

POLYTECH Health & Aesthetics implants are filled with the highly cross-linked EasyFit™ gel, which stands out due to its natural consistency. Because of the gel's "memory effect", it always returns to its original shape after moderate mechanical deformation.

The chemically and mechanically resistant, elastic shell of the implants consists of several layers of a silicone elastomer. In order to maintain a high level of safety, the

shell of the silicone gel-filled implants is equipped with an additional diffusion barrier layer.

The surface of the implant shell is POLYtxt® textured. This surface is known for its reliability and has been proven to contribute to a reduction of capsular contracture rates¹.

Production of the implants is based on a quality management according to ISO 9001 and ISO 13485. All materials used are of implantable medical quality.

1: textured implants contracture rates: 8–15% (Iwuagwu + Frame, 1997), which is 15% less than with smooth implants (Handel, 2006)

Indications

- Breast reconstruction after mastectomy or ablatio mammae
- Contour augmentation and correction
- Adjustment of asymmetries, aplasia or hypoplasia
- Combined breast and thoracic deformities
- Ptosis or atrophy
- Implant exchange

Important Information

The implants are delivered individually packed and sterile.

Attention: Please note the detailed instructions for use included with each product.

POLYTECH Health & Aesthetics GmbH

Altheimer Str. 32 • 64807 Dieburg • Germany

☎ +49 (0) 60 71 98 63 0 • 📠 +49 (0) 60 71 98 63 30

✉ info@polytechhealth.com • www.polytech-health-aesthetics.com

