

# SurgiMend® PRS Meshed

The PRS Family of  
Acellular Dermal Matrices

## Strength and flexibility you can rely on

SurgiMend® PRS Meshed is unique; it is an expandable acellular dermal tissue matrix derived from bovine dermis, which provides up to twice the initial coverage once fully hydrated. Sterile and ready to use ensuring critical and costly theatre time is optimised.

**The first MHRA approved ADM solution for pre-pectoral  
implant based reconstruction.**

**Strong and flexible enough to meet the needs  
of plastic and reconstructive surgeons.**






## What Does It Do?

SurgiMend® PRS Meshed is specifically indicated for plastic and reconstructive surgery, including breast reconstruction. The matrix is used to support and stabilize the breast implant, ensuring it is held in the correct position on the chest wall. It minimizes unwanted movement as SurgiMend® PRS Meshed conforms to the surface of the breast implant more completely than fenestrated products. SurgiMend® PRS Meshed is made from exactly the same bovine dermis we produce our sub pectoral PRS products from and works in harmony with your patients tissues providing sound support – without inflammation.


### SurgiMend® PRS Meshed



Pre-Hydrated Meshed

### Key Advantages

- The ideal solution to pre-pec implant based reconstruction
- High in type III collagen
- Hydrates quickly
- No detrimental inflammatory response
- Classified as a category 4 non irritant by the World Health Organisation
- Excellent handling characteristics with the strength to control and support the implant
- Mesh expansion up to ~18cm x 22cm\*
- Rapid revascularisation to support tissue building
- Easily conforms to round or anatomical breast implants



Post-Hydrated Meshed

### Order Details SurgiMend® PRS Meshed

Reference	Description	Nominal Thickness
SM-606-907-001	20 cm x 10 cm Meshed <i>Expands to ~18 cm x 22 cm*</i>	1 mm

\* When hydrated and fully expanded.

For further information or to place an order, please contact:



#### Q Medical Technologies

Summerlands Trading Estate  
Endmoor, Kendal, Cumbria, LA8 0FB  
T: 0845 194 9284 F: 0845 094 0380  
E: info@qmedical.co.uk W: www.qmedical.co.uk

SurgiMend, Integra and the Integra logo are registered trademarks of Integra LifeSciences in the United States and/or other countries.  
©2017 Integra LifeSciences Corporation. All rights reserved. 0781380-2-EN  
DOCUMENT FOR USE IN EUROPE, MIDDLE-EAST & AFRICA ONLY.

#### Manufacturer :

TEI Biosciences Inc  
7 Elkins Street ■ Boston ■ MA ■ 02127 ■ USA  
Tel: (866) 524-0022 ■ Fax: (888) 623-2259



4-9-11-17V3