


Hands-Free Vein Visualisation, Simplified

Enhance precision, improve workflow, and give clinicians the freedom to focus on what matters most - the patient.



Designed to support faster, more efficient venous access with improved patient comfort, AccuVein's range of hands-free accessories transforms vein visualisation from a handheld tool into a seamless clinical solution.

Whether you're working chairside, bedside, or in a fast-paced environment, our stands and mounting options are designed to support efficiency, consistency, and better patient experiences.



HF580 Powered Stand



HF510 Pole Mount Arm



Accessories

AccuVein offers hands-free solutions and site assessment tools to enhance vein identification, selection and optimal insertion site for successful IV access; also intended to support aesthetic procedures that require vein avoidance.

HF580 Powered Stand

The HF580 is a fully integrated, powered stand designed for hands-free, mobile use across clinical environments.

The HF580 offers:

- ◆ A durable, seamless and flexible arm design that allows for adjustable height positioning, smooth arm rotation and maneuverability to help the clinician make better informed venous site selections and successful needle/catheter insertion.
- ◆ A new and improved break away power cord to improve clinician safety and allow for easy access and storage without worries of tripping or tangling.
- ◆ Durable, ergonomic stand wheels that lock in place by foot activation, allowing the clinician to be truly hands-free to focus on the patient and not the stand.
- ◆ A “quick release” device-holder to easily secure and remove the AV500.

Weight

27 lbs (12 kg), base: 16 lbs (7 kg), arm: 11 lbs (5 kg)

Dimensions

Extended arm reach length: 31" (79 cm)

Base: 22" W x 5" H (56 cm x 13 cm)

Folded height (floor to top of cup): 47" (119 cm)

Carry basket: 13" L x 6" H x 8" W





HF580 Powered Stand: Specifications

Why this matters in practice:

Built-in Charging Capability

No need to use desktop charging cradle

Seamless, Durable Arm Design

Smooth, easy arm rotation for improved manoeuvrability

Quickly Convert to Hands-Free

Allows a single clinician to perform procedures single-handedly

Break-away Power Cord

Safe, easy access and storage
No worries of tangling or tripping

"Quick-Release" Cup Holder

Fast in and fast out, securely held

Carry Basket

Convenient way to transport supplies to patient location

Extended Reach

Arm can reach across patient bed or chair

Stable Base

Compact; easy to move around

Locking Wheels

Foot activated; easy to lock in place

Self-Serviceable

No need to return for repairs

Adjustable Height


Pole can be raised up to 8 inches (20.3 cm)





HF510 Pole Mount Arm

The AccuVein HF510 is a portable, hands-free solution for the AV500 Vein Visualization System.



For the caregiver who doesn't stay in one place, this configuration comes with a clamp that can attach to a rounded bedrail, an IV pole or other flat surfaces; just attach it and your hands are free to perform the procedure.

- A portable arm mount design that is durable, seamless and flexible allowing for adjustable height positioning, smooth arm rotation and maneuverability to help the clinician make better informed venous site selections and successful needle/catheter insertion.
- A "quick release" device-holder to easily secure and remove the AV500.

Key Benefits:

Hands-Free Operation

Perform procedures with complete freedom and control.

Built-In Charging

No need for separate charging stations — always ready when you are.

Extended Reach Arm

Easily position the device across beds or chairs for optimal access.

Smooth, Flexible Movement

Seamless arm rotation for precise positioning.

Stable, Mobile Base

Locking wheels provide both mobility and security.

Integrated Storage Basket

Keep essential supplies within easy reach.



Which Set-Up Is Right For You?

HF580 Powered Stand

- ◆ Fully mobile stand
- ◆ Built-in charging
- ◆ Best for high-volume clinics

HF510 Pole Mount Arm

- ◆ Clamp-based solution
- ◆ No charging
- ◆ Best for flexible / small spaces

Upgrade Your Workflow with Hands-Free Vein Visualisation

From fully integrated powered stands to compact mounting solutions, AccuVein accessories are designed to support clinicians in delivering efficient, precise, and patient-focused care.

Book a demo or speak to our team to find the right setup for your clinic.



Distributed By:

Q Medical Technologies Ltd.
Summerlands Trading Estate
Endmoor, Kendal, Cumbria, LA8 0FB

T: 0845 194 9284

F: 0845 094 0380

E: info@qmedical.co.uk

W: www.qmedical.co.uk